## MIX-25P mini mix instrument



**Product Features** 

Product Numbers : MIX-25P

**Product description**: MIX-30S/MIX-25P/MIX-25 MINI Mixer is used for vibration, blend and stir of sample tissue, cell, bacteria Liquid and chemical reagent in life science and physical and chemical analysis field. Multiple tubes are suitable for 0.2~50ml micro tubes and tubes or small container whose diameter is less than 108mm.

- 1. Mini and compact appearance, powerful function, sucker type under-chassis super shockproof, ideal for high speed operation.
- 2. Multiple mixing modes, touch-operation, continuous operation and timing set.
- 3. Wide speed range:0-3000rpm/min, continuously variable motor and Brushless DC Motor are optional. Stable low speed, powerful high speed.
- 4. Shaking block convenient for replacement, stable and reliable, sturdy and durable eccentric bearing design.
- 5. Various shaking block adapters optional.(Suitable for Eppendorf Tube)

Model	MIX-25 Rotation Type	MIX-25P Rotation Type	MIX-30S Digital display
Operation	Knob + Scale	Knob + Scale	Key + digital display
display pattern			
Orbit	4mm	3 mm	2 mm
Shaking mode	circle	circle	circle
Operation mode	Continuous running/	Continuous running/	Continuous running/
	Short Mix	Short Mix	Short Mix
Speed Range	0 ~ 2500rpm	0 ~ 2500rpm	300 ~ 3000rpm
Timing Range	/	1	1s ~ 999min
Voltage	AC220V/AC110V, 50/60Hz	AC220V/AC110V,	DC24V
		50/60Hz	
Power	60W	60W	30W
Optional Accessories	standard head and flat pad	All	All
Fuse	250V,1A, <b>Φ</b> 5x20	250V,1A, <b>Φ</b> 5x20	1
Dimension	W.132×D.156×H.150mm	W.132×D.156×H.150mm	W.132×D.156×H.150m
			m
Net Weight	2.4kgs	2.4kgs	2.2kgs

## Technical parameters

## Exchangeable blocks for Mix series

Cat#	Description
MX-A	Standard head. Suit for tubes and containers with dia.<30mm
MX-P	Standard platform to hold foam tube holder
MX-B	Convex pad. Use with MX-P. Suit for tubes and containers with
	dia.<99mm
MX-C	Foam tube holder. For 0.2ml tubes, 70 holes. Use with MX-P
MX-D	Foam tube holder. For 0.5ml tubes, 36 holes. Use with MX-P
MX-E	Foam tube holder. For 1.5/2.0ml tubes, 20 holes. Use with MX-P
MX-F	Foam tube holder. For 15ml tubes, 8 holes. Use with MX-P
MX-G	Foam tube holder. For 50ml tubes, 4 holes. Use with MX-P
MX-H	Tube holder platform. For 1.5ml/2.0ml tubes x 22
MX-I	Tube holder platform. For 15ml tubes x 4
MX-J	Tube holder platform. For 50ml tubes x 2
MX-K	Plate holder platform. For all standard plates
MX-L	Universal platform for tube stands MIX-A/B/C and all standard plates
MX-M	Customized

